



## Praxair's *StarWatch Connect*™ Remote Monitoring System

Manage your cryogenic liquid supply from virtually anywhere.


The StarWatch Connect system lets you check the status of all your liquid nitrogen, argon, and carbon dioxide dewars anytime from your computer or mobile device. Users receive highly accurate measurements of dewar liquid level and pressure along with programmable alert notifications by text, email, and phone call.

## **Benefits**

- Help avoid interruptions in cryogenic liquid supply
- Maximize use of product
- Reduce labour needed to check status of dewars
- Better tracking of liquid consumption and plan for dewar replacement







- **1.** StarWatch system provides highly accurate measurements of liquid level and pressure within the cryogenic vessel
  - Liquid level ± 2%
  - Pressure ±1 psig for liquid, ± 5 psig for gas
- 2. Transmitter wirelessly sends liquid level and pressure measurements to Gateway
- 3. Gateway collects measurement data from multiple Transmitters and uploads data to the cloud
- 4. View measurements and alerts on your mobile device
- 5. View measurements and alerts, configure monitoring system, and run reports via website

Optimize cryogenic liquid supply for your critical applications. Contact Praxair today to learn more about the *StarWatch Connect* system.



800.225.8247



specialty\_gas\_support@praxair.com

